



Stability indicating UHPLC method for the assay of maraviroc in bulk and in Formulations

M. Geetha*,

G. Lalitha,

Dr. Shakil Sait

Analytical Research and Development, IPDO, Dr. Reddy's Laboratories Ltd. Bachupally, Hyderabad, Andhra Pradesh, India - 500072.

Department of chemistry, J.N.T.University, Kukatpally, Hyderabad-500072, A.P., India

Corresponding Authors:

Mrs. M. Geetha

Dr Reddys Laboratories Limited, IPDO-Innovation Plaza, Survey No's: 42, 45, 46 & 54, Bachupally, Qutubullapur, Hyderabad, Andhra Pradesh, India - 500 072
E-mail: geetharm@drreddys.com

Abstract:

An Ultra High Pressure Liquid Chromatographic method was developed for the estimation of Maraviroc in bulk and in formulations. The separation was achieved on X bridge (C18 20 x 4.6 mm, 2.5 μ column) using 0.01M potassium dihydrogen phosphate (KH₂PO₄) (pH 7.0 adjusted with ortho phosphoric acid) and acetonitrile (60:40) as mobile phase. The flow rate kept at 0.5 mL/min, column temperature 30°C, and the column eluents were monitored at 210 nm. The forced degradation studies were done to show stability indicating power of the method. The method has been validated accordance with ICH guidelines for specificity, precision, accuracy, linearity, limit of detection, limit of quantification, robustness and ruggedness. The results were found to be well within the limits. The method can be used for the routine analysis of Maraviroc bulk and in formulations.

Keywords: Ultra High Pressure Liquid Chromatography, Maraviroc ICH guidelines, Validation

INTRODUCTION

Maraviroc (UK-427,857) is an orally-available selective chemokine receptor CCR5 antagonist with potent anti-HIV-1 activity. It has exhibited activity against all CCR5-tropic HIV-1 viruses tested, including 43 primary isolates from various diverse sources. Maraviroc was also active against 200 clinically derived HIV-1 envelope-recombinant pseudoviruses, 100 of which were derived from viruses resistant to existing classes.^[1] Maraviroc prevents the binding of the viral envelope, gp120, to CCR5 and, thus, prevents the subsequent membrane fusion events that are necessary for viral entry into the CD41 cells.^[2-3] Maraviroc is not active against CCR2-, CXCR4-and dual-tropic viruses, nor is it cytotoxic.³ Maraviroc binds to the host cells rather than the virus envelope; thus, its mechanism of action differs from enfuvirtide, a

fusion inhibitor.^[4] Additive or synergistic activity has been observed when maraviroc has been assessed in combination with other antiretroviral agents, including abacavir, amprenavir, atazanavir, delavirdine, didanosine, efavirenz, emtricitabine, enfuvirtide, indinavir, lamivudine, lopinavir, nelfinavir, nevirapine, ritonavir, saquinavir, stavudine, tenofovir, zalcitabine, and zidovudine.^[1]

Maraviroc, Chemically 4,4-difluoro-N-[(1S)-3-[(1R,5S)-3-[3-methyl-5-(propan-2-yl)-4H-1,2,4-triazol-4-yl]-8-azabicyclo[3.2.1]octan-8-yl]-1-phenylpropyl]cyclohexane-1-carboxamide is a white to yellowish or brownish powder with a molecular formula of C₂₉H₄₁F₂N₅O and a molecular weight of 513.67. Maraviroc is practically insoluble in water, slightly soluble in ethanol, soluble in methanol, dimethyl sulfoxide and PEG 400.^[4-5]

There are only few liquid chromatographic^[6-8] available for the determination of maraviroc in bulk and formulations. The aim of the present work is to develop a stability indicating cost effective method for estimation of maraviroc in bulk and in formulations.

MATERIALS AND METHODS

Reagents and Chemicals

HPLC gradient grade acetonitrile and methanol from Merck (Mumbai, India) has been used. Potassium dihydrogen phosphate (AR grade), ortho phosphoric acid and triethylamine Solution from Merck have been used. Demineralized water was further purified in the laboratory by filtering through an ultrapure Milli-Q (Millipore, Milford, MA, USA), Maraviroc working standard from Dr. Reddys Laboratories.

Instrumentation and liquid chromatographic conditions

LC was carried out on a Agilent 1290 series UHPLC. The output signal was monitored and processed using chemstation software. The chromatographic column used Waters X bridge C₁₈ 20 x 4.6 mm, 2.5 μ

Standard preparation

Accurately weighed and transferred 100 mg of working standard in to a 100 mL volumetric flask, dissolved and diluted volume with diluent (50:50 water, methanol), filtered through 0.22 μ m filter and sonicated for 10 min. Further 5 ml of solution diluted to 25mL in a 25 mL volumetric flask.

Sample preparation

Accurately weighed and transferred 100 mg of sample in to a 100 ml volumetric flask, dissolved and diluted volume with diluent (50:50 water, methanol), filtered through 0.22 μ m filter and

sonicated for 10 min. Further 5 ml of solution diluted to 25mL in a 25 mL volumetric flask.

Mobile phase preparation

To 1000mL of 0.01M KH₂PO₄, 1mL triethylamine added, pH adjusted to 7.0.with H₃PO₄.

Mixed 600mL of above buffer and 400mL of acetonitrile (v/v) and sonicated for 10 mins.

Optimized Chromatographic Conditions

Column	: X Bridge C ₁₈ 20 x 4.6 mm, 2.5 μ
Column temperature	: 30°C
Wavelength	: 210 nm
Flow rate	: 0.5 mL/min
Injection volume	: 5 μ l
Diluent	: Water and methanol in the ratio of 50:50 v/v

System Suitability

Standard solution was injected six times and % relative standard deviation (%RSD) for peak area (PA) and retention times (RT) calculated. Average of tailing factor (T.F) for individual peak and average theoretical plates (T.P) were calculated. Resolution (RS) between peaks was evaluated.

Specificity

To demonstrate the specificity of the method blank, standard solution and sample solutions were injected.

Degradation studies

To check the performance of the optimized LC method for the separation of degradation products, the drug was subjected to various stress conditions.

Acidic hydrolysis

Maraviroc sample was refluxed with 5N HCl at 60°C for 1hour and then neutralized by adjusting pH to 7.0 with 5N NaOH. The Solution was further diluted to required concentration with diluent.

Note: Sample was not degraded in 0.1N HCL and 1N HCL. So high stress conditions 5N HCL was used.

Alkaline hydrolysis

Maraviroc sample was refluxed with 2N NaOH at 60°C for 1hour and then neutralized by adjusting pH to 7.0 with 2N HCl. The Solution was further diluted to required concentration with diluent.

Note: Sample was not degraded in 0.1N NaOH and 1N NaOH. So high stress conditions 2N NaOH was used.

Oxidative stress

Maraviroc sample was refluxed with 30% 10% H₂O₂ by heating on water bath at 60°C for 1hour. The solution was further diluted to required concentration with diluent.

Photolytic stress

Maraviroc sample was exposed to UV (200watt hour/m²) and Visible (1.2million Lux hours) as per ICH Guidelines. The Sample was prepared as per sample preparation.

Thermal stress

Maraviroc sample was exposed to Temperatures at 105°C for 3days. The Sample was prepared as per sample preparation.

For Water Degradation

Maraviroc sample was refluxed with water by heating on water bath at 100°C for 1hour. The Sample was prepared as per sample preparation.

For Humidity Degradation

Maraviroc sample was exposed to 85% Humidity (Prepared potassium nitrate saturated solution) at 3days. The Sample was prepared as per sample preparation.

Precision

Prepared the standard solution, sample solution prepared at same level at in six replicates, injected into system, chromatograms were recorded, assay of sample solutions was performed and %RSD was calculated.

Intermediate Precision (Ruggedness)

Ruggedness of the method was evaluated by determining the precision of method by analyzing same sample on different system, by different analyst and on different column.

Accuracy

Prepared in triplicate standard solutions at different concentration levels (50-150%), and amount added was calculated in terms of % recovery.

Linearity

Linearity was demonstrated by injecting standard solutions of 10 - 150% with respect to the specification level. Plotted the calibration curve by taking concentration on X-axis and peak area on Y-axis, calculated the correlation coefficient and % y-intercept.

Robustness

Robustness of the method was demonstrated by doing the small variations in mobile phase flow rate, column temperature, organic variation in mobile phase and buffer pH.

RESULTS AND DISCUSSION

Method development

The goal of the method is to detect the impurities and quantify Maraviroc present in bulk drug and in formulations. Based on the pKa value of the drug the pH of the buffer was selected as 7.0. KH₂PO₄ (0.01M), (1mL triethylamine added, pH adjusted to 7.0.with H₃PO₄) was used as buffer. Initially methanol was used as organic phase; the maraviroc peak shape was not good and eluted very lately. Acetonitrile was used as organic phase, initially in 80: 20 buffer: acetonitrile used, peak tailing was observed, acetonitrile content increased to 40% the peak symmetry was good. Zorbax SB C₈ (4.6x50 mm, 1.8μ), Kinetex XB-C-18

(20 x 2.1 mm, 2.5 μ), BEH C₁₈ (20X2.1mm, 2.5 μ) and X Bridge C₁₈ (20 x 4.6 mm, 2.5 μ) columns were tried, good separation with good peak symmetries was observed with X bridge C₁₈ (20 x 4.6 mm, 2.5 μ). At 30°C the peaks symmetry is good and eluted in time. Flow rate was varied between 0.3 mL/min – 0.5 mL/min and 0.5 mL/min flow rate found to be suitable. L injection volumes were tried, at 5 μ L the response was adequate at 5 μ L with good peak shape and size.. On observation 3D spectra of sample in diode array detector response at 210 nm, it has shown good response, hence 210 nm optimized as the detection wavelength.

Method Validation

System suitability

%RSD of retention times, peak areas were less than 1 for maraviroc, average of tailing factor found to be less than 1. Theoretical plates were found to be

more than 5000 for, hence method passes system suitability tests. The results were shown in table no 1.

Table 1: Results of system suitability

Parameter	Result
% RSD of RT	0.32
% RSD of PA	0.39
T. Factor	0.89
T. Plated	9690

Specificity

Blank and Placebo interference

Chromatograms of placebo solutions showed no peaks at the retention time of Maraviroc. This indicates that the excipients used in the formulation do not interfere in the estimation of Maraviroc. The standard and sample chromatograms were identical, peak purity angle, peak purity threshold were good that proves methods specificity.




Fig 1: Blank Chromatogram




Fig 2: Standard Chromatogram




Fig 3: Placebo chromatogram




Fig 4: Sample chromatogram

Interference from degradation products

A study was conducted to demonstrate the effective separation of degradants from Maraviroc and its related known impurities. Separate portions of drug product, drug substance and placebo were exposed to the stress conditions to induce degradation. Stressed samples were injected into the HPLC system with photo diode array detector as per following test method conditions. All degradant peaks were resolved from Maraviroc peak and its related

impurities in the chromatograms of all samples. The Chromatograms of the stressed samples were evaluated for peak purity of Maraviroc using chemstation software. In all forced degradation samples, peak purity was passed for Maraviroc Purity factor is within the purity threshold. The method can be used for determining Maraviroc bulk and pharmaceutical formulations. In all the conditions Maraviroc peak purity angle is less than the purity threshold. Results are shown in table no 2.

S. No	Stress condition	%Net degradation	Peak Purity
1	2 hrs in Water @60°C	4.0	999
2	2 hrs in 2N HCl @ 60°C	3.2	999
3	2 hrs in 3N NaOH @ 60°C	15	999
4	2 hrs in 10%H ₂ O ₂ @ 25°C	3	999
5	Thermal(105°C for 7 days	Nil	999
6	UV and Visible	Nil	999
7	Humidity	Nil	999

Table 2: Results of degradation studies




Fig 5: Water degraded sample




Fig 6: Acid degraded sample




Fig 7: Base degraded sample




Fig 8: Peroxide degraded sample




Fig 9: Photolytic sample




Fig 10: Thermally degraded sample




Fig 11: Humidity sample

Precision and Intermediate Precision

The average amount total impurities found to be 100.4% and % RSD was less 1, when analysis was performed by second analyst on second system the results were well under limits, method is precise.

S. No	Precision	Int. Precision
Preparation-1	100.0	100.4
Preparation-2	100.4	100.0
Preparation-3	100.5	100.7
Preparation-4	100.2	100.7
Preparation-5	101.1	99.7
Preparation-6	100.0	100.6
Mean	100.4	100.4
Std dev	0.4	0.4
%RSD	0.4	0.4

Table 3: Precision and intermediate precision

Accuracy

The mean % recovery values for Maraviroc were found to be 99 – 101%. That proves methods accuracy.

S. No	50%	100%	150%
1	100.1	99.7	99.4
2	99.7	100.1	99.2
3	99.5	100.1	99.3
4	99.9	99.9	99.4
5	100.0	100.7	99.1
6	99.8	99.6	99.3
% RSD	0.2	0.4	0.1

Table 4: Results of accuracy

Linearity

The linearity study reveals that the method is linear from 10 to 150%. The R^2 values for all compounds found to be more than 0.9998.

Conc. in ppm	Area
20	251.39598
50	694.20972
100	1369.08057
140	1904.83203
200	2673.5813
240	3267.24902
300	4054.84863
Corr	0.999890464
Slope	13.52003152
intercept	2.73773641
Bias	0.102399595

Table 5: Results of Linearity




Fig 12: Linearity Graph of Maraviroc

Robustness

There were no changes in results when slight changes in chromatographic conditions were made, method is robust. Results were shown in table no 6.

pH Variation	6.8	7.0	7.2
%RSD of RT	0.37	0.34	0.43
%RSD of PA	0.45	0.37	0.39
T. Factor	0.95	0.91	0.99
T. Plates	9345	9190	8956
Mobile Phase variation	55:45	60:40	65:35
%RSD of RT	0.65	0.37	0.37
%RSD of PA	0.55	0.42	0.34
T. Factor	1.28	0.91	0.89
T. Plates	8745	9428	9634
Temperature Variation	35°C	30°C	25°C
%RSD of RT	0.38	0.42	0.45
%RSD of PA	0.49	0.51	0.42
T. Factor	0.99	0.91	0.97
T. Plates	9190	9712	8990
Flow variation	0.4 mL/min	0.5 mL/min	0.6 ml/min
%RSD of RT	0.44	0.36	0.51
%RSD of PA	0.47	0.45	0.55
T. Factor	0.94	0.89	1.03
T. Plates	8890	9512	9365

CONCLUSIONS

An UHPLC method was developed for the analysis of Maraviroc. The developed method is able to detect maraviroc within 4 min. The

method is sensitive, specific, linear, precise, accurate, robust and easy to perform. This method can be used for the assay of Maraviroc in bulk drug and in formulations.

REFERENCES

- 1) Dorr P, Westby M, Dobbs S, et al. Maraviroc (UK-427,857), a potent, orally bioavailable, and selective small-molecule inhibitor of chemokine receptor CCR5 with broad-spectrum anti-human immunodeficiency virus type 1 activity. *Antimicrob Agents Chemother.* 2005;49(11):4721-4732.
- 2) Rosario MC, Poland B, Sullivan J, Westby M, van der Ryst E. A pharmacokinetic-pharmacodynamic model to optimize the phase IIa development program of maraviroc. *J Acquir Immune Defic Syndr.* 2006;42(12):183-191.
- 3) Westby M, Mori J, Smith-Burchnell C, et al. Maraviroc (UK-427,857)-resistant HIV-1 variants selected by serial passage, are sensitive to CCR5 antagonists and T-20. *Antiviral Ther.* 2005; 10(suppl): S72.
- 4) www.drugbank.ca (Maraviroc)
- 5) www.rxlist.com (Maraviroc)
- 6) V. Kalyana Chakravarthy and D.Gowri Sankar. Stability indicating HPLC method for determination of maraviroc and its degradants/impurities in bulk and pharmaceutical formulation. *Rasayan J. Chem.* 2012; 5 (1): 90-105.
- 7) L. Satyanarayana, S.V. Naidu, M. Narasimha Rao, C. Ayyanna, and Alok Kumar. The Estimation of Maraviroc in Tablet dosage form by RP-HPLC. *Research Journal of Pharmaceutical Dosage Forms and Technology* 2011; 3(5): 230 -235.

Article History:-----

Date of Submission: 22-06-2013

Date of Acceptance: 02-07-2013

Conflict of Interest: NIL

Source of Support: NONE

SJR SCImago
Journal & Country
Rank

Powered by
SCOPUS™